




**NEW!**

# MEB2

(LIQUID GERMICIDAL CONCENTRATE)



Camera's Note

Poor Copy

*Comments*

### MEB2 - LIQUID GERMICIDAL FABRIC TREATMENT PRODUCT FOR LAUNDRIES

MEB2 is a special treatment to be used as an adjunct to standard laundering routine to minimize cross-contamination.

MEB2 aids in preventing or controlling the formation of ammonia.

STANDARD MEB2 aids in preventing mildew for at least 12 days.

#### ACTIVE INGREDIENTS:

- Dimethyl benzyl ammonium chlorides  
neutral (60% C12, 30% C16, 7% C12, 3% C18).....11%
- dimethyl benzyl ammonium chlorides  
neutral (50% C12, 30% C14, 17% C16, 3% C18).....11%
- di-socetyl-paerovethoxyethyl dimethyl  
benzyl ammonium chlorides monohydrate.....30%
- isopropyl alcohol .....20%
- INERT INGREDIENTS .....28%

DANGER: "KEEP OUT OF REACH OF CHILDREN"

SEE BACK PANEL FOR ADDITIONAL CAUTION

...sources of infection... When used as directed... the laundry from... must be rinsed to remove... more than 62% cotton...

WARNING: KEEP OUT OF REACH OF CHILDREN. CAUSES SEVERE EYE AND SKIN DAMAGE. DO NOT GET IN EYES, ON SKIN OR ON CLOTHING. WEAR GOGGLES OR FACE SHIELD AND RUBBER GLOVES WHEN HANDLING. HARMFUL IF SWALLOWED. AVOID CONTAMINATION OF FOOD.

FIRST AID: IN CASE OF CONTACT, IMMEDIATELY FLUSH EYES OR SKIN WITH PLENTY OF WATER FOR AT LEAST 15 MINUTES. FOR EYES, CALL A PHYSICIAN. REMOVE AND WASH CONTAMINATED CLOTHING BEFORE REUSE. IF SWALLOWED, DO NOT INDUCE VOMITING, DRINK LARGE QUANTITIES OF FLUID AND CALL A PHYSICIAN IMMEDIATELY.

STAFGARD M62 is a special treatment to be used as an adjunct to standard laundering routine to minimize cross-contamination.

STAFGARD M62 has been tested by the A.O.A.C. disc-dilution test that can kill gram-positive and gram-negative organisms within a period of 10 minutes.


STAFGARD M62 provides residual benefits, thereby controlling growth of many bacterial encountered during use.

STAFGARD M62 aids in controlling ammonia.

STAFGARD M62 produces a fabric which inhibits the growth of antibiotic resistant Staphylococci and other organisms in the presence of moisture.

STAFGARD M62 does not cause damage of items to discolor when autoclaved.

USDA REG. NO. 7972-3



1.5 Gallon  
Net Weight 8 Pounds

Camerman's Note  
Poor Copy

ACCEPTED WITH COMMENTS