

S/S80

# S/S80 SELF-SANITIZER

FABRIC BACTERIOSTATIC CONCENTRATE  
FOR COMMERCIAL LAUNDRY USE

SOLELY BY  
**RE-ODA CHEM ENGINEERING COMPANY**  
32 WEST ORANGE STREET  
CHAGRIN FALLS, OHIO

S/S80

### FIRST AID


In case of contact, immediately flush eyes or skin with plenty of water for at least 15 minutes. For eyes call a physician. Remove and wash contaminated clothing before re-use. If swallowed, drink promptly large quantities of milk, egg whites, gelatin solution or if these are not available, drink large quantities of water. Avoid alcohol. Call a physician immediately.

### NOTE TO PHYSICIAN:

Probable mucousal damage may contraindicate use of gastric lavage. Measures against circulatory shock as well as oxygen and measures to support breathing manually or mechanically may be needed. If persistent, convulsions may be controlled by the cautious intravenous injection of a short acting barbiturate drug.

Active Ingredients	
n-Alkyl (60% C14, 30% C16, 5% C12, 5% C18) dimethyl benzyl ammonium chlorides	40%
n-Alkyl (68% C12, 32% C14) dimethyl ethylbenzyl ammonium chlorides	40%
Ethyl Alcohol	20%
	<hr/> 100%

E.P.A. REGISTRATION NO. 9687-7



NET CONTENTS

### DANGER:

KEEP OUT OF REACH OF CHILDREN — CORROSIVE

Causes severe skin and eye damage. Do not get in eyes, skin or on clothing. Wear goggles, face shield and rubber gloves when handling. Keep container closed. Wash thoroughly after using. Harmful if swallowed. Avoid contamination of food. Do not reuse empty drum. Return to drum reconditioner or destroy by perforating or crushing and burying in a safe place. **FLAMMABLE!** KEEP AWAY FROM HEAT OR OPEN FLAME. See antidote statement and other precautions on the side panel.

Imparts a self-sanitized finish to the treated fabric that reduces more than 90% of S. aureus and E. coli contained in contaminated wet liquid that may be absorbed in the fabric while in use.

Bacteriostatic against many organisms on treated fabric.

Inhibits the growth of ammonia producing organisms (B. ammoniagenes) which act on the urine on the diaper to produce ammonia.

Eliminates "Static Cling" out of nylon and other synthetics.