

10088-42
U.S. ENVIRONMENTAL PROTECTION AGENCY
OFFICE OF PESTICIDES PROGRAMS
REGISTRATION DIVISION
WASHINGTON, D.C. 20460

4/6/76
EPA REGISTRATION NUMBER

10088-42

DATE OF ISSUANCE

APR 6 1976

NAME OF PESTICIDE PRODUCT

10% LIQUID SANITIZER DISINFECTANT

NOTICE OF REGISTRATION
(Under the Federal Insecticide, Fungicide,
and Rodenticide Act, as amended)

NAME AND ADDRESS OF REGISTRANT

ATHEA LABORATORIES, INC.
4180 - North First St.
Milwaukee, Wisconsin 53212

NOTE: Changes in labeling or formula differing in substance from that accepted in connection with this registration must be submitted to and accepted by the Registration Division prior to use of the label in commerce. In any correspondence on this product always refer to the above EPA registration number.

On the basis of the information furnished by the registrant, the above named pesticide is hereby registered under the Federal Insecticide, Fungicide, and Rodenticide Act.

A copy of the labeling accepted in connection with this registration is returned herewith.

Registration is in no way to be construed as an endorsement or approval of this product by this Agency. In order to protect health and the environment, the Administrator, on his motion, may at any time suspend or cancel the registration of a pesticide in accordance with the Act. The acceptance of any name in connection with the registration of a product under this Act is not to be construed as giving the registrant a right to exclusive use of the name or to its use if it has been covered by others.

Note that this submission was processed under the 1947 Federal Insecticide, Fungicide, and Rodenticide Act. At such time as re-registration is required or amendments are proposed this product will be subject to the Registration, Re-registration and Classification Procedures published in the Federal Register on July 3, 1975. Refer to Section 162.23 of that document. Refer also to PR Notices 75-1 and 75-4.

Joseph M. Tavano
Product Manager (31)
Disinfectants Br.
Registration Division (WH-567)

WH-567:DIS:JMTavano:yh:rm 219:WSME

CC: PED

☐ Attachment is applicable.

SIGNATURE OF APPROVING OFFICIAL

JMT

DATE

6 April 76

COLD DISINFECTION

For disinfection of previously cleaned surgical instruments, barber tools and dental equipment, simply submerge instruments into solution containing 1½ ounce per gallon of water for 10 minutes.

RESTAURANT AND BAR RINSE

Dishes, glassware, silverware, cooking utensils, wash with soap or synthetic detergent; rinse thoroughly and immerse in a sanitizing solution containing 1 ounce to 4 gallons of water. No Terminal Potable Water Rinse Required.

FOR ROOM FOGGING

DIRECTIONS: To minimize the danger of cross-contamination from environmental surfaces apply solution as a fog as described below, before applying standard routine terminal cleaning and disinfecting practices. 1. Remove all human, animal, and plant life from room. 2. Open closet doors and drawers. 3. Set up 34½ revolving platform in center of room. 4. Mount Jet FOGMASTER TRI-JET MODEL 6208, or other fogging device delivering equivalent spray at 3RPM. 5. Fill sprayer reservoir with 1.5 oz. to 1 gallon of water to produce 1150 ppm active ingredients. 6. Set sprayer mechanism to deliver 1 gallon of solution. 7. Fog for 15 minutes for an average 2100 cu. ft. room. For different room sizes, vary spray time proportionately to ensure complete wetting of exposed surfaces. 8. Wait 2 hours before entering room after treatment. 9. Rinse fogging equipment thoroughly with clear water following use.

DANGER

KEEP OUT OF REACH OF CHILDREN

Corrosive, Causes severe eye and skin damage. Do not get in eyes, on skin or on clothing. Protect eyes and skin when handling concentrate. Harmful if swallowed. Avoid contamination of food.

FIRST AID

In case of contact, immediately flush eyes or skin with plenty of water for at least 15 minutes. For eyes, call a physician. Remove and wash all contaminated clothing before reuse. If swallowed, drink milk, egg whites, gelatin solution, or if these are not available, drink large quantities of water. Call a physician.

NOTE TO PHYSICIAN:

Probable musosal damage may contraindicate the use of gastric lavage. Measures against circulatory shock, respiratory depression and convulsion may be needed.

Rinse empty container with water and discard.

DISINFECTANT

SANITIZER

DEODORIZER

ACTIVE INGREDIENTS:

n-Alkyl (60% C14, 30% C16, 5% C12, 5% C18)
dimethyl benzyl ammonium chlorides 5%

n-Alkyl (68% C12, 32% C14) dimethyl
ethylbenzyl ammonium chlorides 5%

INERT INGREDIENTS: 90%
100%

This Product should be diluted before using.

NET CONTENTS: _____

EPA Reg. No. 10088-42-AA-

EPA Est. No. 10088-WI-1

DANGER

KEEP OUT OF REACH OF CHILDREN

See other precautionary and first aid statements on left side panel.

Sold by:

ATHEA LABORATORIES, INC.
4180 NORTH FIRST STREET
MILWAUKEE, WIS. 53212

USE DIRECTIONS

Household

For sanitizing and disinfecting floors, walls and inanimate hard surfaces. For Schools, Homes, Locker Rooms, Garbage Pails, Sink Tops, Corridors, Classrooms, Offices and Shower Stalls. Apply solution with mop or cloth. Sanitize with 1 ounce to 4 gallons of water. Disinfect with 3 ounces to 5 gallons of water. Kills Staphylococcus aureus and Salmonella choleraesuis. For heavily soiled or contaminated areas, a pre-cleaning step is recommended.

Hospitals

For floors, walls and other hard surfaces. In Nursing Homes, Institutions, as well as Sink Tops, Garbage Pails, Telephones and Restrooms, use at 3½ ounces per 5 gallons of water. Apply with cloth, mop, or sponge. At this level will be effective against Pseudomonas aeruginosa PRD-10. For heavily soiled or contaminated areas, a pre-cleaning step is recommended.


Food Processing Equipment

For sanitization of previously cleaned food processing equipment and food utensils dilute 1 ounce per 4 gallons of water to provide 200 ppm of active quaternaries. At this level, No Potable Water Rinse is required.

Dairies

To sanitize dairy equipment such as tanks, lines, pails, and milk cans, first clean and rinse the equipment thoroughly. Then apply sanitizing solution containing 1 ounce to 4 gallons of water (200 ppm). At this level, No Potable Water Rinse is required. Consult Local Health Board for additional information.

This product when used on environmental, inanimate, hard surfaces at 3½ ounces per 5 gallons of water is effective against Influenza A, Herpes, Simplex, Adenovirus type 2, and Vaccinia Viruses.



2/2