

Description Acid Bowl Cleaner is a disinfectant battion is bowl cleaner designed for hospital use The disinfectant toilet bowl cleaner is especially to-inviated to kill many germs commonly found in pospital, industrial, and household environments Disintectant Acid Bowl Cleaner is designed for use in toilet bowls and urinals.

Disinfectant Acid Bowl Cleaner contains an acid compatible detergent system that removes resistant stains and increases its disinfectant action.

When used as directed, Disinfectant Acid Bowl Cleaner cleans, disinfects and deodorizes in one easy step.

## DANGER

### KEEP OUT OF REACH OF CHILDREN

Corrosive. May be fatal or harmful if swallowed. Do not breathe vapor or fumes. Produces chemical burns, Do not get in eyes, on skin or on clothing.

### FIRST AID

If swallowed, call a physician immediately. Drink a teaspoonful or more of milk of magnesia or a cup of scapy water followed with milk or raw egg whites

Eve Contact. Wash with water for 15 minutes. Get prompt medical attention.

Skin Contact. Wipe off the acid gently, immediately. flood the surface with water, using soap freely, then cover with moist magnesia or baking soda

Rinse empty container with water before discarding ∃**t**..

# DISINFECTANT ACID BOWL CLEANER

Designed for Hospital Use Based on Twin Chain Quaternary Ammonium Compounds and Acid Competible Detergent Systems for Hospitals - Industry - Homes Disinfects - Cleans - Deodorizes - Removes Stains

Kills Germs in Toilet Bowls Including Many Found in Hospitals (eg, Pseudomonas aeruginosa & Staphylococcus aureus)

### Active Ingredients

Octyl decyl dimethyl ammonium chloride Dioctyl dimethyl ammonium chloride Didecyl dimethyl ammonium chloride Alkyl<sup>\*</sup> amino betaine Hydrogen chloride

### Inert Ingredients

\*Alkyl=C<sub>8</sub>, 7%; C<sub>10</sub>, 8%, C<sub>12</sub>, 46%; C<sub>14</sub>, 24%; C<sub>16</sub>, 10%; C<sub>18</sub>, 5% Contains No Phosphate Compounds.

### DANGER **KEEP OUT OF REACH OF CHILDREN**

### SEE LEFT PANEL FOR FIRST AID AND ADDITIONAL PRECAUTIONARY STATEMENTS

EPA Reg No 10088-44 AA EPA Est No. 10088-WI 1

### NET CONTENTS

Sold By

e	1.250%
	0.625%
	0.625%
	1.000%
	8.000%
	88.500%
	100.000%

### **USE DIRECTIONS**

Flush toilet. Add 1 fluid ounce (2 tablespoons) directly to water in toilet bowl (% gallon). Brush to cover all bowl surfaces including under the rim. Flush toilet.

For badly stained bowls, repeat above treatment using 2 to 4 fluid ounces of Disinfectant Acid Bowl Cleaner.

For use in urinals, add 1 fluid ounce (2 tablespoons) to 3/4 gallon of water. Apply by brush to all urinal surfaces. Flush.

At 1 fluid ounce/% gallon water, Disinfectant Acid Bowl Cleaner is effective against Pseudomonas aeruginosa, Salmonella chloeraesuis and Staphylococcus aureus (in toilet bowls and urinals). This performance has been confirmed by the AOAC Use-Dilution Test Procedure.

The cationic detergent system of Disinfectant Acid Bowl Cleaner inhibits the acid corrosion of metal parts found around toilet booths and urinals.

ACCEPTED 1 of an end of a strategy with the first terms of THE REPORT OF POLICE FOR USE A REPAIR NO. STATECT TO FT ACHE COMMENTS الا المتعلمية بمتهميهم مستعمر والروام المتعادة المتعاملة والمار المراجع

\$ 91