


# DYNA-KLEEN™ -A

HEAVY DUTY CLEANING CONCENTRATE GERMICIDE



**USE DIRECTIONS:**

Use at 2 ounces per gallon. This product will be effective against *Pseudomonas aeruginosa* PRD-10 at this level. This has been substantiated by the AOAC Use Dilution Confirmation Test. At this level, the product is fungicidal against pathogenic fungi.

For best results, use Dyna-Kleen-A with a cloth or mop and apply to walls, floors and other hard surfaces. No rinse is required.


If used for sanitization of food equipment or food contact items, limit the level of active quaternary to 200 ppm (2.8 ounces per 5 gallons).

This product when used on environmental, inanimate, hard surfaces at 2 ounces per gallon of water, is effective against Influenza A<sub>2</sub>, Herpes Simplex, Adenovirus type 2, and Vaccinia Viruses.

**PHENOL COEFFICIENTS:**

Guaranteed	Minimum
S. aureus.....	42
S. typhosa.....	27

Rinse empty container thoroughly with water and discard.



CLEANER  
DEODORIZER

DISINFECTANT  
FUNGICIDE

Designed for Hospital and Institutional use.

**ACTIVE INGREDIENTS:**

- n-Alkyl (60% C<sub>14</sub>, 30% C<sub>16</sub>, 5% C<sub>12</sub>, 5% C<sub>18</sub>) dimethyl benzyl ammonium chlorides .....2.25%
- n-Alkyl (68% C<sub>12</sub>, 32% C<sub>14</sub>) dimethyl ethylbenzyl ammonium chlorides .....2.25%
- Sodium Carbonate .....3.00%
- Tetrasodium ethylenediamine tetraacetate ....1.00%

INERT INGREDIENTS: .....91.50%  
100.0%

Contains no phosphorous compounds. The surfactants in this product are biodegradeable.

NET CONTENTS: 1 Qt. 32 Oz. EPA Reg. No. 10612-3

**WARNING**

Keep Out of Reach of Children. See First Aid and other precautions on right panel.

Manufactured by:

**TZL LABORATORIES**  
East Syracuse, New York 13057

Dyna-Kleen-A is a cleaner, disinfectant, and sanitizer, designed for institutional use such as Hospitals, Nursing Homes and Institutions where housekeeping will be of prime importance in controlling the hazard of cross-infection. Dyna-Kleen-A can provide the assurance you need in this area, when used as directed. For use on floors, walls, sink tops, garbage pails, telephones, rest rooms, bed pans and operating rooms. Effective against *Staphylococcus aureus*, *Salmonella choleraesuis* and *Pseudomonas aeruginosa*. Also effective against pathogenic fungi, *Trichophyton interdigitale*.

When used for sanitization of food processing equipment and food contact items, limit active quaternary to 200 ppm. Treated food contact surfaces must be rinsed with potable water prior to reuse.

**WARNING**

Keep Out of Reach of Children. Causes eye and skin irritation. Do not get in eyes, on skin or on clothing. Harmful if swallowed. Avoid contamination of food.

**FIRST AID**

In case of contact, immediately flush eyes or skin with plenty of water for at least 15 minutes. For eyes, call a physician. Remove and wash all contaminated clothing before reuse. If swallowed, drink milk, egg whites, gelatin solution, or if these are not available, drink large quantities of water. Call a physician.