

DIRECTIONS FOR USE

General Classification

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

M133 detergent/disinfectant has been designed specifically for hospitals, nursing homes, schools, and other institutions where housekeeping is of prime importance.

DISINFECTANT-To disinfect hard inanimate surfaces (such as walls, floors, table tops), add 1 ounce M133 detergent/disinfectant per gallon of water. Apply solution with mop, cloth, or sponge or mechanical sprayer so as to wet thoroughly. Allow to remain wet for 10 minutes and then let air dry.

To disinfect toilet bowls, flush toilet, add 1 ounce M133 detergent/disinfectant directly to the bowl water. Swab the bowl completely using a scrub brush or toilet mop, making sure to get under the rim. Let stand for 10 minutes and flush.

For heavily soiled areas, a pre-cleaning step is required and recommended. Prepare a fresh solution for each use.

MILDEWSTAT-To control mold and mildew on pre-cleaned, hard non-porous surfaces (such as floors, walls, table tops) add 1 ounce M133 per gallon of water. Apply solution with a cloth, mop, or sponge making sure to wet all surfaces completely. Let air dry. Prepare a fresh solution for each use. Repeat application at weekly intervals or when mildew growth reappears.

PRECAUTIONARY STATEMENTS

Hazard to Human and Domestic Animals
DANGER


Keep Out of Reach of Children. Causes eye and skin irritation. Do not get in eyes, skin or on clothing. Harmful if swallowed.

STATEMENT OF PRACTICAL TREATMENT

In case of contact, immediately flush eyes or skin with plenty of water for at least 10 minutes. For eyes, call a physician. Remove and wash all contaminated clothing before reuse. If swallowed, drink milk, egg whites, gelatin solution or if there are not available, drink large quantities of water. Call a physician.

NOTE TO PHYSICIAN

Probable mucosal damage may contraindicate the use of gastric lavage. Measures against circulatory shock, respiratory depression and convulsion may be needed.



M133 detergent/disinfectant

heavy duty

hospital institutional *of (M133)*

CLEANER

DETERGENT

DEODORIZER

GERMICIDE

MILDEWSTAT (on hard inanimate surfaces)

DISINFECTANT

FUNGICIDE (against

pathogenic fungi)

VIRUCIDE*

ACTIVE INGREDIENTS:

n-Alkyl (80%C14, 30%C18, 5%C12, 5%C18)	
dimethyl benzyl ammonium chlorides	4.50%
n-Alkyl (80%C12, 3%C14)	
dimethyl ethylbenzyl ammonium chlorides	4.50%
INERT INGREDIENTS	91.00%
Total	100%

EPA REG NO. ~~47772-4~~ 47772-4
EPA EST. NO.

ONLY FOR SALE TO, USE, AND STORAGE BY SERVICE PERSONS

KEEP OUT OF REACH OF CHILDREN

DANGER

see side panels for additional precautionary and first aid statements

MFST
CHEMICAL COMPANY
4802 Birchtree Lane
Temple Hills, MD 20748

M133 is a phosphate free concentrated detergent/disinfectant which provides effective cleaning, deodorizing, disinfection, and ~~disinfection~~. M133 detergent/disinfectant has been designed specifically for hospitals, nursing homes, schools, and other institutions where housekeeping is of prime importance in controlling the hazard of cross contamination. M133 when diluted at the rate of 1 ounce per gallon of water is an effective disinfectant against the organism Pseudomonas aeruginosa and meets all requirements for hospital use. When used as directed, M133 is formulated to disinfect inanimate hard surfaces such as walls, floors, sink tops, toilet bowls, tables, chairs, telephones and bed frames. For larger areas such as operating rooms, patient care facilities and restrooms, the product is designed to provide both general cleaning and disinfecting. In addition, M133 deodorizes those areas which generally are hard to keep fresh smelling, such as garbage storage areas, empty garbage bins and cans, toilet bowls and any other areas which are prone to odors caused by microorganisms. **DISINFECTANT**-At 1 ounce per gallon dilution, M133 exhibits effective disinfectant activity against the following:

Pseudomonas aeruginosa PRD-10

Escherichia coli	Brevibacterium ammoniogenes
Klebsiella pneumoniae	Streptococcus faecalis
Salmonella schultzei	Shigella dysenteriae
Salmonella choleraesuis	Enterobacter aerogenes
Staphylococcus aureus	Staphylococcus aureus

At this level the product is also fungicidal against the pathogenic fungi, Trichosporon mentagrophytes

VIRUCIDAL* ACTIVITY: The product when used on environmental inanimate hard surfaces at 1 ounce per gallon of water exhibits effective virucidal* activity against Influenza A2/Asian (representative of the common flu virus), Herpes simplex (causative agent of fever blisters and mononucleosis), Adeno virus type 2 (causative agent of upper respiratory infections), and Vaccinia virus (representative of the pox virus group).

EFFICACY TESTS HAVE DEMONSTRATED THAT M133 IS AN EFFECTIVE BACTERICIDE, FUNGICIDE, AND VIRUCIDE* IN THE PRESENCE OF ORGANIC SOIL (5% BLOOD SERUM), MOLD and MILDEWSTAT. M133 prevents and controls mold and mildew and the odors they cause on hard non-porous inanimate surfaces.

STORAGE AND DISPOSAL

Pesticide: Discard any excess pesticide by flushing down drain with large quantity of water.
Container: Do not reuse empty container. Rinse empty container with water and put in trash collection.

BEST AVAILABLE COPY

NET CONTENTS

Reg # 47772-4
P.M. 31

P. 1/1