

BEST DOCUMENT AVAILABLE

America's  
Finest  
Quality  
Filtration  
Systems



**Aqua Guard Bacteriostatic  
Under-Counter Water Filter  
Replacement Cartridge  
Model U-602R**

ACTIVE INGREDIENT SILVER (as elemental) 00.026%  
INERT INGREDIENT ACTIVATED CARBON 99.974%

**CAUTION: KEEP OUT OF REACH OF CHILDREN.**

**DISPOSAL:** Wrap spent filter cartridge in newspaper and  
discard with trash

**CAPACITY:** 600 gallons

**FLOW RATE:** One-half gallon per minute set manually

**CONTENTS:** One Bacteriostatic Water Filter

**SEE ENCLOSED CIRCULAR FOR ADDITIONAL INFORMATION**

EPA EST. No. APPLIED FOR

EPA REG. No. APPLIED FOR

AQUA GUARD, Inc., Rt. 15, Box 588, Tyler, Texas 75707



**CAUTION**  
Unit must be installed  
with water flowing the  
same direction as  
arrows

ACCEPTED  
with conditions  
by EPA L...

JUL 23 1988

U.S. DEPARTMENT OF HEALTH AND HUMAN SERVICES  
FEDERAL BUREAU OF INVESTIGATION  
RECEIVED  
EPA REG. No. 50402-5

## AQUA GUARD BACTERIOSTATIC UNDERCOUNTER WATER FILTER REPLACEMENT CARTRIDGE Model U-602R

The Aqua Guard Bacteriostatic Undercounter Water Filter Replacement Cartridge Model U-602R is uniquely designed to provide you with clearer, good-tasting and good-smelling drinking water. It inhibits the growth of bacteria entrapped within the filter medium. The "effective-life-expectancy" of this Replacement Cartridge is 600 gallons; or 12 months of normal use. Replacement of the cartridge in the Aqua Guard Bacteriostatic Undercounter Water Filter Model U-602 is very simple if you will follow these directions.

**DIRECTIONS FOR USE:** It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

1. Turn on cold water and open valve on filtered water faucet allowing water to run into sink. Slowly turn the "T" bar on the self-tapping valve (attached to cold water line under counter-see Fig. "A") until water has stopped flowing from filtered water faucet into sink. Double-check to be sure water will not flow through filtered faucet and has been completely shut off at self-tapping valve.

2. Release the spent cartridge from the cartridge holder bracket.

3. Holding the spent cartridge in one hand, turn the nut on the plastic compression fitting located on top of the unit (see Fig. "B") counterclockwise until red vinyl tubing may be easily removed. **NOTE:** place pan under sink to catch excess water when tubing is disconnected.

4. Repeat same procedure with compression nut on bottom of spent cartridge, releasing blue vinyl tube (see Fig. "C"). Wrap spent cartridge in newspaper and discard with trash.

5. Loosen compression nut on top of new Replacement Cartridge. Insert red vinyl tube into compression fitting and tighten by hand only.

6. Repeat same procedure with compression nut on bottom of new Replacement Cartridge. Insert blue vinyl tube, and tighten by hand only.

7. Place new Replacement Cartridge into cartridge holder bracket and secure.




8. With valve on filtered faucet still open and cold water turned on, slowly open self-tapping valve until unit fills and water runs out filtered faucet into sink. Slowly close self-tapping valve, measuring flow with a container, until flow-rate of 1/2 gallon-per-minute is reached. **Important** for maximum quality filtration.

9. Let filter run for approximately 10 minutes, watching under counter for leaks. (Do not drink this water). Slightly tighten any fitting that leaks.

10. Complete warranty card and mail immediately.

This filter cartridge is designed for use with the Aqua Guard Bacteriostatic Undercounter Water Filter Model U-602. It should only be used in conjunction with municipally-treated tap water in the home.

**AQUA GUARD, INC. • Rt. 15, Box 588 • Tyler, Texas 75707 • (214) 561-4545**



BEST DOCUMENT AVAILABLE

JUL 28 1984

57400

<b>U.S. ENVIRONMENTAL PROTECTION AGENCY</b> <b>OFFICE OF PESTICIDES PROGRAMS</b> REGISTRATION DIVISION (WH-567) WASHINGTON, D.C. 20460	EPA REGISTRATION NO.	DATE OF ISSUANCE 23 JUL 1984
	TERM OF ISSUANCE	
	NAME OF PESTICIDE PRODUCT	
<b>NOTICE OF PESTICIDE:</b> <input type="checkbox"/> REGISTRATION <input checked="" type="checkbox"/> REREGISTRATION <i>(Under the Federal Insecticide, Fungicide, and Rodenticide Act, as amended)</i>		
NAME AND ADDRESS OF REGISTRANT (Include ZIP code)		
BEST DOCUMENT AVAILABLE		
<b>NOTE:</b> Changes in labeling formula differing in substance from that accepted in connection with this registration must be submitted to and accepted by the Registration Division prior to use of the label in commerce. In any correspondence on this product always refer to the above U.S. EPA registration number.		
On the basis of information furnished by the registrant, the above named pesticide is hereby Registered/Reregistered under the Federal Insecticide, Fungicide, and Rodenticide Act.		
A copy of the labeling accepted in connection with this Registration/Reregistration is returned herewith.		
Registration is in no way to be construed as an indorsement or approval of this product by this Agency. In order to protect health and the environment, the Administrator, on his motion, may at any time suspend or cancel the registration of a pesticide in accordance with the Act. The acceptance of any name in connection with the registration of a product under this Act is not to be construed as giving the registrant a right to exclusive use of the name or to its use if it has been covered by others.		
<input type="checkbox"/> ATTACHMENT IS APPLICABLE		
SIGNATURE OF APPROVING OFFICIAL		DATE

U.S. ENVIRONMENTAL PROTECTION AGENCY  
OFFICE OF PESTICIDES PROGRAMS  
REGISTRATION DIVISION (77-567)  
WASHINGTON, D.C. 20460

EPA REGISTRATION NO.

DATE OF ISSUANCE

50402-5

23 JUL 1984

TERM OF ISSUANCE

NOTICE OF PESTICIDE:  REGISTRATION  
 REREGISTRATION  
(Under the Federal Insecticide, Fungicide,  
and Rodenticide Act, as amended)

NAME OF PESTICIDE PRODUCT

Aqua Guard Bacteriostatic Under-  
counter Water Filter Replacement  
Cartridge Model U-602R

NAME AND ADDRESS OF REGISTRANT (Include ZIP code)

Aqua Guard, Inc.  
Rt. 15, Box 588  
Tyler, TX 75707

NOTE: Changes in labeling formula differing in substance from that accepted in connection with this registration must be submitted to and accepted by the Registration Division prior to use of the label in commerce. In any correspondence on this product always refer to the above U.S. EPA registration number.

On the basis of information furnished by the registrant, the above named pesticide is hereby Registered/Reregistered under the Federal Insecticide, Fungicide, and Rodenticide Act.

A copy of the labeling accepted in connection with this Registration/Reregistration is returned herewith.

Registration is in no way to be construed as an indorsement or approval of this product by this Agency. In order to protect health and the environment, the Administrator, on his motion, may at any time suspend or cancel the registration of a pesticide in accordance with the Act. The acceptance of any name in connection with the registration of a product under this Act is not to be construed as giving the registrant a right to exclusive use of the name or to its use if it has been covered by others.

This product is conditionally registered in accordance with FIFRA sec.

3(c)(7)(A) provided that you:

1. Submit and/or cite all data required for registration/reregistration of your product under FIFRA sec. 3(c)(5) when the Agency requires all registrants of similar products to submit such data.

2. Make the labeling changes listed below before you release the product for shipment:

a. Add the phrase "EPA Registration No. 50402-5."


b. Include the EPA Establishment Number on the label where it is referenced.

The silver release data methods of chemical analysis for silver are acceptable. For your records, please note that the data have been filed under EPA Accession Number 253933.

3. Submit five (5) copies of your final printed labeling before you release the product for shipment. Refer to the A-79 Enclosure for a further description of final printed labeling.

If these conditions are not complied with, the registration will be subject to cancellation in accordance with FIFRA sec. 6(e). Your release for shipment of the product constitutes acceptance of these conditions.

A stamped copy of the label is enclosed for your records.



John E. Lee  
Product Manager (31)  
Disinfectants Branch  
Registration Division (TS-767C)

DISCOUNTS APPLICABLE

SIGNATURE OF APPROVING OFFICIAL

DATE

10F4  
Pm 21  
- 20-5

America's  
Finest  
Quality  
Filtration  
Systems



**Aqua Guard Bacteriostatic  
Under-Counter Water Filter  
Replacement Cartridge  
Model U-602R**

ACTIVE INGREDIENT: SILVER (as elemental) 00.026%  
INERT INGREDIENT: ACTIVATED CARBON 99.974%

**CAUTION: KEEP OUT OF REACH OF CHILDREN.**

**DISPOSAL:** Wrap spent filter cartridge in newspaper and discard with trash.

**CAPACITY:** 600 gallons.

**FLOW RATE:** One-half gallon per minute, set manually.

**CONTENTS:** One Bacteriostatic Water Filter.

**SEE ENCLOSED CIRCULAR FOR ADDITIONAL INFORMATION**

EPA EST. No. APPLIED FOR

EPA REG. No. APPLIED FOR

AQUA GUARD, Inc., R.L. 15, Box 588, Tyler, Texas 75707

ACCEPTED  
with COMMENTS  
in EPA Letter Dated:

JUL 23 1984

Under the Federal Insecticide,  
Fungicide, and Rodenticide Act  
as amended, for the pesticide  
registered under EPA Reg. No.

50402-5



water filter model U-602 is very simple if you will follow these directions.

**DIRECTIONS FOR USE:** It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

1. Turn on cold water and open valve on filtered water faucet allowing water to run into sink. Slowly turn the "T" bar on the self-tapping valve (attached to cold water line under counter-see Fig. "A") until water has stopped flowing from filtered water faucet into sink. Double-check to be sure water will not flow through filtered faucet and has been completely shut off at self-tapping valve.

2. Release the spent cartridge from the cartridge holder bracket.

3. Holding the spent cartridge in one hand, turn the nut on the plastic compression fitting located on top of the unit (see Fig. "B") counterclockwise until red vinyl tubing may be easily removed. NOTE: place pan under sink to catch excess water when tubing is disconnected.

4. Repeat same procedure with compression nut on bottom of spent cartridge, releasing blue vinyl tube (see Fig. "C"). Wrap spent cartridge in newspaper and discard with trash.




5. Loosen compression nut on top of new Replacement Cartridge. Insert red vinyl tube into compression fitting and tighten by hand only.

6. Repeat same procedure with compression nut on bottom of new Replacement Cartridge. Insert blue vinyl tube, and tighten by hand only.

7. Place new Replacement Cartridge into cartridge holder bracket and secure.

8. With valve on filtered faucet still open and cold water turned on, slowly open self-tapping valve until unit fills and water runs out filtered faucet into sink. Slowly close self-tapping valve, measuring flow with a container, until flow-rate of 1/2 gallon-per-minute is reached. Important for maximum quality filtration.

9. Let filter run for approximately 10 - 15 minutes.



BEST DOCUMENT AVAILABLE